## THERMO DISTANCE

Termometro Infrarossi Frontale a distanza

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
## EN Infrared No Contact Forehead Thermometer THERMO DISTANCE

### Legend:

- LCD display
- Proximity light
- Power (ON/OFF) and measurement (SCAN) button
- Operating mode (MODE) button
- Memory (MEM) button
- Infrared sensor
- Proximity sensor
- Battery compartment (2 x 1.5V "AAA")

## WARNINGS AND USEFUL TIPS

⚠ Attention – In any circumstances, the reading provided by this thermometer must only be considered as a reference. Any medical treatments and/or interventions must only be taken after consulting a doctor.

The Chicco Thermo Distance thermometer is designed for the intermittent measurement of the human body's temperature in a domestic environment, and can be used upon people of all ages.

The following recommendations must be carefully observed during the product's use. Any activities that are inconsistent with or do not take into account these recommendations could result in personal injury and/or could affect the accuracy of the thermometer itself.

- ⚠ Attention – Unauthorized modifications to the device are not permitted. Do not modify the product in any way without the manufacturer's prior authorization.
- In order to avoid incorrect measurements, make sure that the lens on the thermometer's infrared sensor is clean and intact prior to use.
- Avoid touching the lens directly with your fingers
- It is recommended to perform at least three consecutive temperature measurements. If the readings should differ from one another, the average value should be considered as the reference temperature.
- Do not expose the thermometer to extreme temperatures or humidity levels: follow the instructions provided in this manual. Do not expose to direct sunlight.
- Avoid dropping the device or subjecting it to strong jolts.
- Both the patient and the thermometer should be kept within the same room under stable environmental conditions for at least 30 minutes prior to each measurement.
- Avoid taking temperature measurements for at least 30 minutes after physical activity, bathing,

- swimming, consuming food and/or hot/cold beverages or spending time outdoors.
- Wait at least one minute between one temperature measurement and the next. Measurements performed too frequently could result in incorrect temperature readings. Otherwise, use the average value of the readings as a reference.
  - There is no absolute body temperature value. Keep reliable records of your own body temperature as a reference for determining whether or not you have a fever.
  - Contact your local distributor for any assistance you may require with configuration, use and maintenance operations.
  - The batteries must be kept strictly out of the reach of children, as their ingestion could result in poisoning and/or other serious health risks.
  - Remove all traces of sweat, make-up, creams, etc. from the forehead a few minutes prior to performing any temperature measurements.
  - If the infrared sensor is pointed in the direction of a hot object, it will perform the temperature measurement remotely. In order to avoid incorrect readings, the infrared sensor must not be pointed at any heat sources, and the interval that passes from the time in which the thermometer is ready and the time in which the forehead is scanned must be kept to a minimum.
  - If possible, the measurements should always be taken by pointing

the infrared sensor at the centre of the forehead.

- Holding the thermometer in your hands for too long and/or exposing the device to external heat sources could result in distorted temperature readings. For this reason, the body temperature reading could result as being higher or lower than the actual value.
- Any use upon children or individuals with limited cognitive abilities should be supervised by an adult.

## **MEASURING THE TEMPERATURE OF THE FOREHEAD**

Remember that the thermometer needs to be left in the room in which the measurement is to be taken for at least 30 minutes prior to use.

1. Any attempts at taking measurements in areas of the body other than the forehead could provide inaccurate results.
2. Both the child and the adult must remain still during the measurement.
3. Forehead temperature measurements taken using infrared thermometers are normally equivalent to oral temperature measurements taken using digital or liquid metal thermometers.
4. Temperature measurements carried out during sleep should not be directly compared with measurements during wakefulness, because the temperature during sleep is generally lower.

	<p>Press the power button. The thermometer will carry out a self-test and all the symbols will simultaneously appear on the display.</p>
	<p>An acoustic signal will sound and the screen will display the symbols shown in the figure. Make sure that the thermometer is set to forehead measurement mode "👤".</p>
	<p>Hold down the "SCAN" button to perform the measurement. The "measurement in progress" "🔋" symbol will appear on the display.</p>
	<p>Holding down the "SCAN" button, point the infrared sensor to the centre of the forehead. The proper distance position is indicated by the proximity light turning on and by the sounding of a continuous intermittent acoustic signal.</p>
	<p>Release the "SCAN" button. After about one second, an acoustic "beep" will sound, indicating that the temperature has been measured. At this point, the reading will appear on the display.</p>
	<p><b>FEVER CHECK</b>  <math>&lt; 38^{\circ}\text{C}/100.4^{\circ}\text{F}</math> = for temperatures below <math>38^{\circ}\text{C}/100.4^{\circ}\text{F}</math>, the screen will display a smiley face symbol  <math>\geq 38^{\circ}\text{C}/100.4^{\circ}\text{F}</math> = for temperatures above or equal to <math>38^{\circ}\text{C}/100.4^{\circ}\text{F}</math>, the screen will display a sad face symbol, followed by three beeps</p>
	<p>After 30 seconds, the thermometer will emit an acoustic signal and will automatically shut itself off.</p>

5. Clean the thermometer's lens after each use. In order to ensure the accuracy of the reading, wait at least 10 minutes before performing another measurement.

## MEASURING THE TEMPERATURE OF OBJECTS/LIQUIDS

	Press the power button. The thermometer will carry out a self-test and all the symbols will simultaneously appear. The device will then emit two beeps and the display will show the symbols "00" and "👤".
	In order to switch to "object measurement" mode, press the "MODE" button once. The screen will display the symbols shown in the figure.
	Hold down the "SCAN" button to perform the measurement. The "measurement in progress" symbol "◀■■■■" will appear on the display.
	While holding down the "SCAN" button, point the infrared sensor to the object holding the thermometer at approximately 4~6cm away.
	Release the "SCAN" button. After 1 second, the device will emit a single beep and the temperature reading will appear on the display.
	After 30 seconds, the thermometer will emit an acoustic signal and will automatically shut itself off.

## NIGHT FUNCTION

The innovative "NIGHT FUNCTION" allows children temperature to be measured while sleeping, without waking them up.

With the thermometer on, press the "MODE" button to select the desired operating mode

Forehead function    Object function    Forehead/night function    Object/night function

## MEMORY FUNCTION

	With the thermometer on, press the memory button (MEM) twice. The display will show the last measurement performed, accompanied by the "M" symbol.
	The following symbols will appear alongside each memory item: "M" indicating the body temperature "M" with a thermometer icon indicating the object temperature
	Press the ON/MEM button repeatedly to view the memory items in sequence. This thermometer is capable of storing a series of up to 10 of recent readings.
	Hold down the "MEM" button for at least 3 seconds to clear the records stored in the device's memory. An acoustic signal will sound and the screen will display the symbols shown in the figure.

## SELECTING THE MEASUREMENT SCALE (°C/°F)

	The thermometer is configured to perform measurements in Celsius (°C) by default. In order to change the unit of measure from °C to °F, make sure that the thermometer is on and press the "MODE" and "MEM" buttons simultaneously. The selected unit of measure will be saved for any subsequent measurements. The same procedure can be used to change the unit of measure from °F back to °C.
--	--

## BATTERY REPLACEMENT

	When the battery's charge is low, the battery icon will appear on the display. Although the thermometer can continue to be used, the battery will have to be replaced as soon as possible. Once the battery has been entirely drained, the "Lo" symbol will appear on the display as well.
	Open and remove the battery compartment's cover. Remove the drained batteries and insert 2 new 1.5V "AAA" alkaline batteries into the battery compartment, respecting the proper polarity as shown in the figure. Reposition the battery compartment's cover.

## ⚠ WARNINGS FOR REPLACING THE BATTERIES AND PROPER DISPOSAL OF USED BATTERIES

1. Always use new batteries of the type and specifications indicated in this manual: do not use old batteries together with new batteries or with batteries of different types.
2. Do not use rechargeable batteries, as these may be of inferior quality and duration. The use of rechargeable batteries could compromise the device's performance.
3. Dispose of the used batteries promptly, keeping them out of the reach of children and away from any heat sources.
4. For environmental protection purposes, the used batteries must be disposed of at appropriate collection sites and in compliance with the current local or national laws, policies and/or regulations. The batteries utilized comply with the European directive 2006/66/EC, and the relative warnings must be observed for their disposal.
5. The batteries must be kept strictly out of the reach of children, as their ingestion could result in poisoning and/or other serious health risks.
6. It is recommended to remove the batteries whenever the device is not expected to be used for an extended period of time.
7. The device does not have a shelf life and has a typical service life of three years. The batteries used inside the device, however, have a shelf life of 6 months (half a year): electrolyte leakage could take place once the batteries' shelf life has expired. The batteries must be replaced once their foreseen service life has been reached (1000 measurements). If in doubt, or to request clarifications, please contact your local distributor.

## CARE AND MAINTENANCE

	<b>MEASUREMENT LENS/SENSOR</b> Clean with a swab soaked in alcohol. Never clean the lens directly with water.
	<b>THERMOMETER</b> Clean the thermometer with a soft dry cloth. Do not use abrasive products or solvents. Take care to avoid scratching the lens or the display. Do not expose the thermometer to direct sunlight, jolts and/or extreme temperature or humidity levels.

## TROUBLESHOOTING

In the event of a malfunction or incorrect measurement, a relative error message will appear.

Error message	Problem	Solution
	The detected temperature is lower than: 1. Forehead measurement mode: 10.0°C (50.0°F) 2. Object measurement mode: 0°C (32.0°F)	Only use the thermometer within the specified temperature range. Clean the sensor if necessary.
	The detected temperature is greater than: 1. Forehead measurement mode: 50.0°C (122.0°F) 2. Object measurement mode: 100°C (212.0°F)	If the problem persists, contact your nearest retailer.
	The operating temperature is not in the range from 16°C ~ 35°C (60.8°F ~ 95.0°F)	Only use the thermometer within the specified temperature range.

## TECHNICAL FEATURES

- Measuring range [°C]  
10.0°C – 50.0°C (50.0°F – 122.0°F) [Human body]  
0°C – 100°C (32.0°F – 212.0°F) [Object]
- OPERATION DISTANCE 4 ~ 6cm (±1cm)
- Measuring accuracy HUMAN BODY 22°C – 40°C (71.6°F – 104.0°F): ±0.3°C (0.5°F)
- Measuring accuracy OBJECT T ≤20°C: ±1.0°C; T >20°C: ±5%
- Repeatability HUMAN BODY ±0.15 °C (0.3°F) according to ASTM E1965-98.

- Measuring Time HUMAN BODY & OBJECT Around 1 second
- Operating Environment 16°C ~35°C (60.8°F – 95°F) with relative humidity up to 95% (non condensing)
- Storage & Transportation Environment

 -25°C ~ 55°C (-13°F – 131°F)

-  with relative humidity up to 95% (non condensing).

## CONFORMITY STANDARDS

This products complies with the essential requirements of MDD directive 93/42/EE C . The following standards have been applied during product design and manufacturing:

- ASTM E1965-98 Standard Specification for Infrared Thermometers for Intermittent Determination of Patient Temperature

- IEC/EN 60601-1 Medical electrical equipment -- Part 1: General requirements for basic safety and essential performance
- IEC/EN 60601-1-2 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests

## EMC Tables

The Chicco Infrared No Contact Forehead Thermometer REF. 00006931000000 is intended for use in the electromagnetic environment specified below. The customer or the user of the Chicco Infrared No Contact Forehead Thermometer REF. 00006931000000 should assure that it is used in such an environment.

Guidance and manufacturer's declaration-electromagnetic emissions			
Emissions test	Compliance	Electromagnetic environment-guidance	
RF emissions	CISPR 11	The Chicco Infrared No Contact Forehead Thermometer REF. 00006931000000 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF emissions CISPR 11	Class B	The Chicco Infrared No Contact Forehead Thermometer REF. 00006931000000 is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.	
Harmonic emissions IEC 61000-3-2	Not applicable		
Voltage fluctuations Flicker emissions IEC 61000-3-3	Not applicable		
Guidance and manufacturer's declaration-electromagnetic immunity			
Immunity test	IEC 60601 Test level	Compliance Level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Floor should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Radiated RF IEC 61000-4-3	3V/m 80MHz to 2.5 GHz	$E_1 = 3 \text{ V/m}$	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Chicco Infrared No Contact Forehead Thermometer REF. 00006931000000, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance.</p> $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey," should be less than the compliance level in each frequency range." Interference may occur in the vicinity of equipment marked with the following symbol (Ⓜ).</p>
------------------------------	-----------------------------	-----------------------	--

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations from radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast can not be predicted theoretic call with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Chicco Infrared No Contact Forehead Thermometer REF. 00006931000000 is used exceeds the applicable RF compliance level above, the Chicco Infrared No Contact Forehead Thermometer REF. 00006931000000 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Chicco Infrared No Contact Forehead Thermometer REF. 00006931000000.

b. Over the frequency range 150kHz to 80MHz, field strengths should be less than  $[V1] \text{ V/m}$ .

**Recommended separation distances between portable and mobile RF communications equipment and the Chicco Infrared No Contact Forehead Thermometer REF. 00006931000000**

The Chicco Infrared No Contact Forehead Thermometer REF. 00006931000000 is intended for use in an electromagnetic environment in which radiated RF distances are controlled. The customer or the user of the Chicco Infrared No Contact Forehead Thermometer REF. 00006931000000 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Chicco Infrared No Contact Forehead Thermometer REF. 00006931000000 as recommended below according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150kHz to 80MHz $d=1.2 \sqrt{P}$	80kHz to 800MHz $d=1.2 \sqrt{P}$	800kHz to 2.5GHz $d=2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 8- MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**SYMBOLS KEY**

CE marking (Notified body 0197 - TÜV Rheinland LGA Products GmbH)

Refer to instruction

Type BF applied part

Operating instructions

Warning

**IP 22**

Protection level against solid foreign objects and dripping water (15° tilted)

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### **This product complies with the Directive**

**2012/19/EU.** The crossed bin symbol on the appliance indicates that the product, at the end of its life, must be disposed of separately from domestic waste, either by taking it to a separate waste disposal site for electric and electronic appliances or by returning it to your dealer when you buy another similar appliance. The user is responsible for taking the appliance to a special waste disposal site at the end of its life. If the disused appliance is collected correctly as separate waste, it can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and health, and contributes towards the recycling of the product's materials. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

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**THIS PRODUCT COMPLIES WITH EC DIRECTIVE 2006/66/EC.** The crossed bin symbol on the batteries or product pack indicates that,

at the end of their life, they must not be disposed of as urban refuse. They must be disposed of separately from domestic waste, either by taking them to a separate waste disposal site for batteries or by returning them to your dealer when you buy similar rechargeable or non-rechargeable batteries. The chemical symbols

Hg, Cd, Pb, printed under the crossed bin symbol, indicate the type of substance contained in the batteries: Hg=Mercury, Cd=Cadmium, Pb=Lead. The user is responsible for taking the batteries to a special waste disposal site at the end of their life, so that they can be treated and recycled. If the spent batteries are collected correctly as separate waste, they can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and human health, and contributes towards the recycling of the batteries' substances. Noncompliance with the norms on battery disposal damages the environment and human health. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the batteries.

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## WARRANTY

This thermometer is guaranteed free of manufacturing faults, under normal conditions of use according to the manual instructions, for a period of up to two years from the date of purchase. This guarantee does not apply to damage resulting from improper use. Should repair be needed within the warranty period, send the thermometer with the proof of purchase (which must be provided) to your distributor, who will repair it or replace it free of charge. This warranty does not give the purchaser any particular legal rights; the rights of the purchaser may vary from country to country or state to state. In case of problems not reported in this manual, please contact your Art-sana S.p.A dealer or distributor. The product used and maintained according to instructions has an unlimited life span. In the event of breakage or damage, dispose of according to instructions.

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AVITA Corporation
Wellkang Ltd.,
Artsana S.p.A. - Via Saldarini Catelli 1 - 22070 Grandate (CO) Italy